

OCI-BQ-HD180L-06

130.00 mm

60.00 mm

Hemodialyzer • Hämodialysator • Hémodialyseur • Hemodializador • Hemodialisador • Emodializzatore • Hemodialysator • Hemodialysator • 透析器 • 血液透析器 • 血液透析器

**Polyethersulfone Hollow
Fiber Hemodialyzer**
Hollow Fiber Dialyzer

Medi Sense

OCI-HD180L

Low Flux

Membrane Material : PES

Membrane Area : 1.8m²




Chengdu OCI Medical Devices Co., Ltd.
No.2401, West Port Avenue, Southwest Airport Economic
Development Zone, Shuangliu District, Chengdu, Sichuan
Province, 610299, China

Tel: +86-28-67085899 Fax: +86-28-67085880

EU REP Lepu Medical (Europe) Coöperatief U.A.
Abbe Lenstra Boulevard 36, 8448 JB, Heerhooft, The Netherlands

Tel: +31-515-573399 Fax: +31-515-760020
Importer: Superior Healthcare International House # 132 Street # 08
Al Moor Garden Madinat Tawon Faisalah
ELK-01389 Registration No. MDR-2000001448 M.R.P. Rs.3000

Maximum TMP: 500mmHg(66.5kPa)



Use only with an ultrafiltration control machine
Store in a dry place with ventilated, clean,
non-corrosive air within the temperature of 0~40°C,
relative humidity ≤80%.

COM52Y39K0

C82M56Y4K0

请仔细核对文字内容

2026.1.13